



Swiss Waagen DC GmbH

Freiestr. 13, 8610 Uster
 Fax 043 843 95 92
 Telefon 043 843 95 90

www.swisswaagen.ch
 shop@swisswaagen.ch

Kalibriergewichte gem. OIML-Fehlergrenzen

Nennwert	E1 +/-	E2 +/-	F1 +/-	F2 +/-	M1 +/-	M2 +/-	M3 +/-
1 mg	0.003 mg	0.006 mg	0.020 mg	0.060 mg	0.020 mg		
2 mg	0.003 mg	0.006 mg	0.020 mg	0.060 mg	0.020 mg		
5 mg	0.003 mg	0.006 mg	0.020 mg	0.060 mg	0.020 mg		
10 mg	0.003 mg	0.008 mg	0.025 mg	0.080 mg	0.025 mg		
20 mg	0.003 mg	0.010 mg	0.03 mg	0.10 mg	0.3 mg		
50 mg	0.004 mg	0.012 mg	0.04 mg	0.12 mg	0.4 mg		
100 mg	0.005 mg	0.016 mg	0.05 mg	0.16 mg	0.5 mg	1.6 mg	
200 mg	0.006 mg	0.020 mg	0.06 mg	0.20 mg	0.6 mg	2.0 mg	
500 mg	0.008 mg	0.025 mg	0.08 mg	0.25 mg	0.8 mg	2.5 mg	
1 g	0.010 mg	0.03 mg	0.10 mg	0.3 mg	1.0 mg	3.0 mg	10 mg
2 g	0.012 mg	0.04 mg	0.12 mg	0.4 mg	1.2 mg	4.0 mg	12 mg
5 g	0.016 mg	0.05 mg	0.16 mg	0.5 mg	1.6 mg	5.0 mg	16 mg
10 g	0.020 mg	0.06 mg	0.20 mg	0.6 mg	2.0 mg	6.0 mg	20 mg
20 g	0.025 mg	0.08 mg	0.25 mg	0.8 mg	2.5 mg	8.0 mg	25 mg
50 g	0.03 mg	0.10 mg	0.3 mg	1.0 mg	3.0 mg	10 mg	30 mg
100 g	0.05 mg	0.16 mg	0.5 mg	1.6 mg	5.0 mg	16 mg	50 mg
200 g	0.10 mg	0.3 mg	1.0 mg	3.0 mg	10 mg	30 mg	100 mg
500 g	0.25 mg	0.8 mg	2.5 mg	8.0 mg	25 mg	80 mg	250 mg
1 kg	0.5 mg	1.6 mg	5.0 mg	16 mg	50 mg	160 mg	500 mg
2 kg	1.0 mg	3.0 mg	10 mg	30 mg	100 mg	300 mg	1 000 mg
5 kg	2.5 mg	8.0 mg	25 mg	80 mg	250 mg	800 mg	2 500 mg
10 kg	5.0 mg	16 mg	50 mg	160 mg	500 mg	1 600 mg	5 000 mg
20 kg	10 mg	30 mg	100 mg	300 mg	1 000 mg	3 000 mg	10 g
50 kg	25 mg	80 mg	250 mg	800 mg	2 500 mg	8 000 mg	25 g
100 kg		160 mg	500 mg	1 600 mg	5 000 mg	16 g	50 g
200 kg		300 mg	1 000 mg	3 000 mg	10 g	30 g	100 g
500 kg		800 mg	2 500 mg	8 000 mg	25 g	80 g	250 g
1000 kg		1 600 mg	5 000 mg	16 g	50 g	160 g	500 g